

- Drug free
- Non-invasive
- Achieve an erection in minutes¹

Drug-free
Success
that works



Rapport™

For more information please visit us at owenmumford.com

Our accreditations include:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



1. Yuan, J., Hoang, A., Romero, C., Lin, H., Dai, Y. and Wang, R. (2010). Vacuum therapy in erectile dysfunction—science and clinical evidence. International Journal of Impotence Research, 22(4), pp.211-219.

© Copyright Owen Mumford 2017
RAPPORT/GLUTIMCC/0917/01

Rapport™

Rapport™

Drug-free treatment for erectile dysfunction

Rapport™ vacuum therapy devices have been developed as user-friendly, effective treatments for erectile dysfunction. The advantages of using Rapport include:

- A non-drug therapy, which can be used with other medications.¹
- Applied externally and therefore non-invasive.
- An erection can be achieved within minutes and therefore spontaneity for intercourse is not compromised.¹

Rapport™ Classic

Rapport™ Classic is a unique vacuum therapy system developed as an effective treatment for erectile problems. This vacuum therapy system creates a vacuum around the penis which quickly produces a satisfactory erection by drawing blood into the penis. The pump is then removed leaving the constriction ring at the base of the penis.

The Rapport™ Classic system consists of:

- 1 clear plastic vacuum tube
- 1 hand pump with a pressure relief valve
- 5 sizes of penile rings
- 2 sizes of transfer sleeve(s)
- 1 loading cone
- A tube of water-soluble lubricating jelly
- A full set of patient instructions
- All the above items are supplied in a discreet carrying case



Rapport™ Rings

The **Rapport™ Rings** loading system is ideal for men who need help to maintain an erection. The constriction ring, placed at the base of the penis, will maintain an erection until it is removed.



Rapport™ Rings system consists of:

- 2 sizes of transfer sleeve(s)
- 1 loading cone
- 5 sizes of penile rings
- A tube of water-soluble lubricating jelly
- A full set of patient instructions