










Unilet® Excelite

		International	North American
	Unilet® Excelite 23G	Box of 100: 001155	Box of 100: AT0545
	Unilet® Excelite 26G	Box of 100: 001153	Box of 100: AT0555
	Unilet® Excelite 28G	Box of 100: AT0535 Box of 100: AT0935 (UK only)	Box of 100: AT0565
	Unilet® Excelite 30G	Box of 100: 000405	Box of 100: AT0575
	Unilet® Excelite 33G	Box of 100: AT9585	Box of 100: AT0585

Unilet® ComforTouch™

	Unilet® ComforTouch™ 23G		Box of 100: AT0455
	Unilet® ComforTouch™ 26G		Box of 100: AT0435
	Unilet® ComforTouch™ 28G	Box of 100: AT0465 Box of 200: AT0460	Box of 100: AT0925
	Unilet® ComforTouch™ 30G		Box of 100: AT0465

For more information please visit us at owenmumford.com

Our accreditations include:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



1. Simbec Research. A single-blind, randomised, 8-way crossover study to compare the blood volume and pain perception of capillary blood sampling. UK; 2007 p. 1-15.
 2. Owen Mumford. User Evaluation Report 15-011. 2015. 3. Except Roche Multiclix and Softclix Plus, LifeScan OneTouch Delica, Nipro TrueDraw. 4. Owen Mumford
 User Evaluation Report 14/002 and 14/004.2014. 5. Owen Mumford. Unilet Excelite LCR. 2014. 6. Owen Mumford (2015) LCR DWR000127. UK

© Copyright Owen Mumford 2017
 AT/AUTO-UNI-RNG-HCP-LT/MCC/08/17/04

Flexible finger sampling



Autolet®
Lancing devices

Unilet®
Single use lancets

Autolet® impression

Lancing device

Autolet® Impression is a comprehensive lancing device designed to provide maximum comfort and control when obtaining a blood sample.



Precise adjustment for a tailored experience.

To achieve optimal blood flow and comfort, **Autolet® Impression** offers 21 settings, allowing you to tailor the experience. Seven depth settings and three force settings allow exact adjustment for all skin types and sampling requirements.

Maximum Comfort.

The unique raised dots of **Comfort Zone Technology®** mask the sensation of pain, making for a comfortable lancing experience¹.




Maximum Comfort

Comfort Zone Technology® minimises the pain of finger sampling¹. When pressed against your finger, a series of eight raised dots on the platform stimulate the nerve endings. This sends a powerful message of comfort to the brain, masking the weaker pain stimulus from the needle.



Masking the
pain
sensation



Autolet® Plus

Advanced lancing device

Autolet® Plus is an intuitive lancing device which has been designed to ensure obtaining a blood sample is comfortable and easy.



Flexible.

With **Autolet® Plus** you can choose from nine depth settings to be confident in obtaining sufficient blood for glucose testing.

Easy to use.

Autolet® Plus automatically primes during lancet insertion and has a simple mechanism for release, making the device easy to use.

Unilet®


Single use lancets

Two ranges of single use lancets are available, **Unilet® Excelite** and **Unilet® ComforTouch™**. The range of lancet designs and needle specifications allows you to select the most suitable lancet for your needs. Both lancets are compatible with most capillary blood lancing devices^{3,5,6}.

Unilet® Excelite

Comfort.

The **Unilet® Excelite** lancet provides a comfortable lancing experience without compromising on blood volume⁴.



Compatibility.

The **Unilet® Excelite** lancet will fit most capillary blood lancing devices making it simple to choose a lancet to meet your needs^{3,5}.

Unilet® ComforTouch™

Triple sharpened.

The **Unilet® ComforTouch™** lancets have triple sharpened steel needle tips designed for maximum comfort.



Convenient.

You can push the needle into the needle cap for convenient storage before disposal.

100% needle inspection guarantees superior quality.

Every **Unilet® ComforTouch™** lancet is inspected to ensure the highest possible standards of quality.

Compatibility.

The **Unilet® ComforTouch™** lancet will fit most capillary blood lancing devices^{3,6}.