

HEAD OFFICE

Owen Mumford Ltd
Brook Hill, Woodstock
Oxfordshire OX20 1TU, United Kingdom
T : +44 (0)1993 812021
E : info@owenmumford.co.uk

UK / IRELAND SALES & MARKETING

Owen Mumford Ltd
Primsdown Industrial Estate
Worcester Road, Chipping Norton
Oxfordshire OX7 5XP, United Kingdom
T : +44 (0)1608 645555
E : info@owenmumford.co.uk

NORTH AMERICA

T : +1 770-977-2226
E : info@owenmumfordinc.com

FRANCE

T : +33 (0)2 32 51 88 70
E : fr.info@owenmumford.com

GERMANY

T : +49 (0)6026/97750
E : de.info@owenmumford.com

CHINA

T : +86 (0)21 6032 3528
E : cn.info@owenmumford.com

MALAYSIA

T : +60 (0)7 553 9617
E : my.info@owenmumford.com

MEXICO

T : +52 (55) 3733 8759
E : mx.info@owenmumford.com

owenmumford.com

OUR ACCREDITATIONS INCLUDE:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



Controlled self injection delivering for you

 **OWEN MUMFORD**

Autoject[®] 2

Self-injection devices



More reason to smile

owenmumford.com

Autoject® 2

Self-injection devices

Many people are apprehensive about injecting themselves using syringes. **Autoject® 2** has been developed to make it a more manageable and simple experience. Psychological barriers are overcome by hiding the syringe and needle away, whilst the whole process is made simple through improved ergonomics and one-handed operation.

One-handed operation increases choice of injection sites.

Autoject® 2 makes it possible to use a wider range of injection sites than manually injecting with a syringe alone, e.g. when injecting in difficult to reach sites such as into the back.

Easy-to-use design makes it possible for more people to self-inject.

The ergonomic, easy-to-use design allows people, especially children and seniors, to self-inject.

Safety interlock mechanism.

When the **Autoject® 2** is loaded with a syringe, the safety device prevents accidental firing until the device is pressed against the skin.

Hidden needle reduces anxiety.

Autoject® 2 conceals the needle from the user's view prior to injecting, reducing anxiety for patients who are frightened by the sight of needles.

Other key features

- Visible indicator to show when injection is complete.
- Safety mechanism to prevent accidental firing.
- Depth adjustment to tailor injection.
- Compatible with a wide range of syringe brands.

Autoject® 2 is not for use with glass syringes.

Autoject® 2 Range

Autoject® 2 Fixed needle

Compatible with most major brands of fixed needle syringes.



Autoject® 2 Removable needle

Specifically designed for compatibility with the 1ml BD Tuberculin syringe.



Autoject® 2 automates the process of self-injection. At the touch of a button **Autoject® 2** injects the needle to a pre-set depth and delivers the contents of the syringe.

There are two standard **Autoject® 2** models:

1. For patients using fixed needle syringes.
2. For patients using removable needle syringes.

Each **Autoject® 2** comes in a zipped soft pouch that includes:

- 1 **Autoject® 2**
- 2 depth adjusters (fixed needle)
1 depth adjuster (removable needle)
- a yellow primer/depth guide
- instructions for use



Fixed needle Autoject® 2 set shown above