

Available by Prescription



 **Amielle®** Comfort  
Graduated Vaginal Dilators

## Help Overcome Vaginismus and Dyspareunia



Amielle® Comfort is a vaginal dilator kit designed to help alleviate vaginal discomfort and penetration problems caused by vaginismus and dyspareunia<sup>1</sup>, including discomfort following childbirth and gynecological surgery.

*Amielle dilators are available by prescription and (in many cases) are covered by private insurance.*

### Amielle Comfort Vaginal Dilator System

CPT Code: Procedure 57400 & 58999 | Supplies Code: 99070 | HCPCS Code: E1399

AMCOMSELL/OMI/0417/1/US 1. Data on file.

## About Amielle Comfort

Amielle Comfort vaginal dilators are designed to assist women in the treatment of vaginal conditions such as vaginismus and dyspareunia<sup>1</sup>, including discomfort following childbirth, surgery or radiotherapy.

## Vaginismus & Dyspareunia

These vaginal dilators can be an effective aid to complement sex therapy treatment of vaginismus, and help overcome painful sexual intercourse associated with dyspareunia.

## Following Surgery or Childbirth

Amielle Comfort dilators can also be used to gradually reduce pain or anxiety sometimes experienced during intercourse following surgery or childbirth. Offered in five graduated sizes for a natural progression of treatment, the largest dilator is designed to represent the size of an average erect penis to help patients resume normal sexual activity.



## PRODUCT DETAILS



### Dilators Also Available Separately



Size 1: 2.75"x.6"  
SM2101 | NDC 08470-2101-01



Size 2: 3.54"x.8"  
SM2102 | NDC 08470-2102-01



Size 3: 4.33"x1"  
SM2103 | NDC 08470-2103-01



Size 4: 5.51"x1.2"  
SM2104 | NDC 08470-2104-01



Size 5: 6.30"x1.4"  
SM2105 | NDC 08470-2105-01

### Rx Only Vaginal Dilator Kit Includes:

- 5 graduated exerciser cones
- Universal handle
- Lubricating gel
- Easy-to-read instruction guide
- Discreet carrying bag

SM2100 | HCPCS: E1399 | NDC 08470-2100-01

CPT Code: Procedure 57400 & 58999 | Supplies Code: 99070 | HCPCS Code: E1399

AMCOMSELL/OM/0417/1/US 1. Data on file.



For additional information,  
please contact Customer Service at 1-800-421-6936  
or at [info@owenmumfordinc.com](mailto:info@owenmumfordinc.com)