

**HEAD OFFICE**

Owen Mumford Ltd  
Brook Hill, Woodstock  
Oxfordshire OX20 1TU, United Kingdom  
T : +44 (0)1993 812021  
E : info@owenmumford.co.uk

**UK / IRELAND SALES & MARKETING**

Owen Mumford Ltd  
Primsdown Industrial Estate  
Worcester Road, Chipping Norton  
Oxfordshire OX7 5XP, United Kingdom  
T : +44 (0)1608 645555  
E : info@owenmumford.co.uk

**NORTH AMERICA**

T : +1 770-977-2226  
E : info@owenmumfordinc.com

**FRANCE**

T : +33 (0)2 32 51 88 70  
E : fr.info@owenmumford.com

**GERMANY**

T : +49 (0)6026/97750  
E : de.info@owenmumford.com

**CHINA**

T : +86 (0)21 6032 3528  
E : cn.info@owenmumford.com

**MALAYSIA**

T : +60 (0)7 553 9617  
E : my.info@owenmumford.com

**MEXICO**

T : +52 (55) 3733 8759  
E : mx.info@owenmumford.com

[owenmumford.com](http://owenmumford.com)

**OUR ACCREDITATIONS INCLUDE:**

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



© Copyright Owen Mumford 2016  
AH/OP/6000-LTMCC/0116/02



Developed to make eye  
drop application easy

 OWEN MUMFORD

**AutoDrop**<sup>®</sup>

Eye drop guide



More reason to smile

**AutoSqueeze**<sup>™</sup>

Bottle squeezer



[owenmumford.com](http://owenmumford.com)

# AutoDrop<sup>®</sup>



Eye drop guide

Taking eye drops can be a challenge for many patients, young or old. **AutoDrop<sup>®</sup>** and **AutoSqueeze<sup>™</sup>** have been developed to make self administration as easy as possible.

**AutoDrop<sup>®</sup>** and **AutoSqueeze<sup>™</sup>** are an eye drop treatment reusable for the length of the treatment, making eye drop application simple.

## Unique design prevents blinking.

A pinhole draws the patient's gaze away from the descending drop, whilst a small lip holds the lower eyelid down, to prevent blinking.



## Delivers the drop to the right place.

The nozzle is positioned perfectly over the eye ensuring that drops are successfully delivered.



Developed together with the Royal National Institute for the Blind.

# AutoSqueeze<sup>™</sup>



Combine **AutoSqueeze<sup>™</sup>** with **AutoDrop<sup>®</sup>**, to make squeezing of the bottle easier.

Many patients who use eye drops may have other conditions that affect their dexterity and ability to squeeze small eye drop bottles. **AutoSqueeze<sup>™</sup>** clips onto the bottle to provide a secure, ergonomic grip and extra leverage. It can be used in conjunction with **AutoDrop<sup>®</sup>** or on its own.

## Proven compatibility with most eye drop bottles.

**AutoDrop<sup>®</sup>** has been tested with a wide range of eye drops for compatibility<sup>1</sup>, for a full, up-to-date list please contact Owen Mumford.

1. Owen Mumford (2010) Autodrop Compatibility Testing, NPD\_REP\_0413 Issue: 1. UK.

**AutoDrop<sup>®</sup>** and **AutoSqueeze<sup>™</sup>** have been designed to overcome common obstacles to patient compliance:

- Missing the eye with the drop.
- Blinking as the drop is delivered.
- Difficulty in squeezing the bottle to deliver the drop.